



Human IL23A recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5123
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Pro189
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	18.7kDa
FORM	
BUFFER	PBS,pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9NPF7
ALIAS/SYNONYMS	SGRF

