



## Human PCSK9 recombinant protein

### OVERVIEW

---

|                           |                            |
|---------------------------|----------------------------|
| <b>PRODUCT REFERENCE</b>  | PX-P5124                   |
| <b>ORIGIN SPECIES</b>     | Homo sapiens (Human)       |
| <b>FRAGMENT</b>           | Gln31-Gln152&Ser153-Gln692 |
| <b>TAG</b>                | C-terminal His Tag         |
| <b>DELIVERY CONDITION</b> | Dry Ice                    |

### PRODUCT INFORMATION

---

|                                |  |
|--------------------------------|--|
| <b>EXPRESSION SYSTEM</b>       | Eukaryotic expression  |
| <b>HOST</b>                    | Mammalian cells  |
| <b>MOLECULAR WEIGHT</b>        | 72.82kDa   |
| <b>FORM</b>                    |  |
| <b>BUFFER</b>                  | PBS,pH7.5  |
| <b>STABILITY &amp; STORAGE</b> | 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection) |
| <b>UNIPROT</b>                 | Q8NBP7   |
| <b>ALIAS/SYNONYMS</b>          | NARC1  |

