



Human CD319 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5125
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Ser23-Ser225
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	22.33kDa
FORM	
BUFFER	PBS,pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9NQ25
ALIAS/SYNONYMS	CS1, SLAMF7,CD2-like receptor-activating cytotoxic cells, CRACC, Membrane protein FOAP-12, Novel Ly9,Protein 19A

