



Human CD72 (Arg117-Asp359) recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5126
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Arg117-Asp359
TAG	Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	26.73kDa
FORM	
BUFFER	PBS,pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P21854
ALIAS/SYNONYMS	

