



CSPG4 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5182
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Ser1705-Ile1942
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	27.34 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q6UVK1
ALIAS/SYNONYMS	chondroitin sulfate proteoglycan 4, Chondroitin sulfate proteoglycan NG2, MCSP

