



CCL11 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5186
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Gly24-Pro97
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	8.14 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P51671
ALIAS/SYNONYMS	SCYA11, eotaxin-1, eotaxin1, eosinophil chemotactic protein

