



HBV-C S recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5192
ORIGIN SPECIES	Hepatitis B virus genotype C (isolate Vietnam/3270/2000)
FRAGMENT	Met175-Ile400
TAG	
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	25.32 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9E6S4
ALIAS/SYNONYMS	Large envelope protein, L glycoprotein, L-HBsAg, LHB, Large S protein, Large surface protein, Major surface antigen, Hepatitis B virus

