



STX2 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5194
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Arg188
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	24.14 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P32856
ALIAS/SYNONYMS	STX2C, STX2B, STX2A, Syntaxin-2, Epimorphin, STX2, EPIM

