



## CD62L recombinant protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P5195
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	Glu109-Ser346
<b>TAG</b>	N-GST Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	Escherichia coli (E.coli)
<b>MOLECULAR WEIGHT</b>	52.82 kDa
<b>FORM</b>	
<b>BUFFER</b>	PBS, pH7.5
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	P14151
<b>ALIAS/SYNONYMS</b>	CD62 antigen-like family member L, Leukocyte adhesion molecule 1, CD62L, LYAM1, TQ1, SELL, Leukocyte-endothelial cell adhesion molecule 1, L-selectin, LAM-1, gp90-MEL, LNHR, LECAM1, Lymph node homing receptor, Leukocyte surface antigen Leu-8

