



HIV1 gp120 / SU recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5202
ORIGIN SPECIES	HIV 1
FRAGMENT	Ala29-Arg513
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	57.06 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9YUL5
ALIAS/SYNONYMS	Surface protein gp120, Glycoprotein 120

