



TNFRSF12A recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5208
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Glu28-Trp79
TAG	C-terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	30.72 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9NP84
ALIAS/SYNONYMS	FN14, TWEAKR

