



Bakers yeast GPM1 / PGAM 1, C-His, recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5453
ORIGIN SPECIES	Bakers yeast
FRAGMENT	Met1-Lys246
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	28.58 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P00950
ALIAS/SYNONYMS	

