



ATP5IF1, N-GST, recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5507
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Asp106
TAG	N-terminal GST Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	39.09 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9UII2
ALIAS/SYNONYMS	

