



CCL26 / Eotaxin-3 / TSC-1, N-GST, recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5540
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Ser27-Leu94
TAG	N-terminal GST Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	34.91 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9Y258
ALIAS/SYNONYMS	

