



CCL3 / MIP-1-alpha, N-GST, recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5543
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Ala27-Ala92
TAG	N-terminal GST Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	34.29 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P10147
ALIAS/SYNONYMS	

