



CD158b2 / KIR2DL3, C-Fc, recombinant protein

OVERVIEW

| | |
|---------------------------|----------------------|
| PRODUCT REFERENCE | PX-P5574 |
| ORIGIN SPECIES | Homo sapiens (Human) |
| FRAGMENT | Met1-His245 |
| TAG | C-terminal Fc Tag |
| DELIVERY CONDITION | Dry Ice |

PRODUCT INFORMATION

| | |
|--------------------------------|--|
| EXPRESSION SYSTEM | Eukaryotic expression |
| HOST | Mammalian cells |
| MOLECULAR WEIGHT | 53.27 kDa |
| FORM | |
| BUFFER | PBS, pH7.5 |
| STABILITY & STORAGE | 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection) |
| UNIPROT | P43628 |
| ALIAS/SYNONYMS | |

