



CD262 / TNFRSF10B / TRAIL-R2, C-Fc, recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P5602
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Glu182
TAG	C-terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	46.35 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O14763
ALIAS/SYNONYMS	

