



CXCL16 / SR-PSOX, N-GST, recombinant protein

OVERVIEW

| | |
|---------------------------|----------------------|
| PRODUCT REFERENCE | PX-P5679 |
| ORIGIN SPECIES | Homo sapiens (Human) |
| FRAGMENT | Pro49-Pro137 |
| TAG | N-terminal GST Tag |
| DELIVERY CONDITION | Dry Ice |

PRODUCT INFORMATION

| | |
|--------------------------------|--|
| EXPRESSION SYSTEM | Prokaryotic expression |
| HOST | Escherichia coli (E.coli) |
| MOLECULAR WEIGHT | 36.93 kDa |
| FORM | |
| BUFFER | PBS, pH7.5 |
| STABILITY & STORAGE | 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection) |
| UNIPROT | Q9H2A7 |
| ALIAS/SYNONYMS | |

