



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

FGL2, N-His, recombinant protein

OVERVIEW

| | |
|---------------------------|----------------------|
| PRODUCT REFERENCE | PX-P5720 |
| ORIGIN SPECIES | Homo sapiens (Human) |
| FRAGMENT | Ile209-Pro439 |
| TAG | N-terminal His Tag |
| DELIVERY CONDITION | Dry Ice |

PRODUCT INFORMATION

| | |
|--------------------------------|--|
| EXPRESSION SYSTEM | Prokaryotic expression |
| HOST | Escherichia coli (E.coli) |
| MOLECULAR WEIGHT | 29.27 kDa |
| FORM | |
| BUFFER | PBS, pH7.5 |
| STABILITY & STORAGE | 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection) |
| UNIPROT | Q14314 |
| ALIAS/SYNONYMS | |



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.