



Cellular tumor antigen p53

OVERVIEW

PRODUCT REFERENCE	PX-P6000
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Asp393
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	46.96kDa
FORM	
BUFFER	TBS, pH 8.0, 0.02% Sarcosyl
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P04637
ALIAS/SYNONYMS	Cellular tumor antigen p53, Tumor suppressor p53, Tp53

