



Recombinant Monkeypox virus/MPXV B2R/Hemagglutinin Protein, N-His

OVERVIEW

PRODUCT REFERENCE	PX-P6009
ORIGIN SPECIES	Monkeypox virus (MPXV)
FRAGMENT	Ser19-Asp274
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	30.23 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

UNIPROT

ALIAS/SYNONYMS

