



Recombinant Monkeypox virus/MPXV H3L Protein, N-His

OVERVIEW

PRODUCT REFERENCE	PX-P6011
ORIGIN SPECIES	Monkeypox virus (strain Zaire-96-I-16) (MPX)
FRAGMENT	Lys77-Leu279
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	25.86 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	
ALIAS/SYNONYMS	

