



Human PD-L2 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6014
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Thr220
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	23.62 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9BQ51
ALIAS/SYNONYMS	PDCD1 ligand 2, B7-DC, PD-L2, PDCD1LG2, PDCD1L2, PDL2, Butyrophilin B7-DC, Programmed death ligand 2, Programmed cell death 1 ligand 2, CD273, B7DC, PD-1 ligand 2

