



Human CD113 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6019
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Thr171-Ser356
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	22.82 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9NQS3
ALIAS/SYNONYMS	CDw113, PRR3, PVRL3, Poliovirus receptor-related protein 3, NECTIN3, Nectin cell adhesion molecule 3, CD113, Nectin-3

