



Human CD209 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6020
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Glu70-Ala404
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	40.61 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9NNX6
ALIAS/SYNONYMS	CD209, DC-SIGN1, Dendritic cell-specific ICAM-3-grabbing non-integrin 1, CD209 antigen, CLEC4L, DC-SIGN, C-type lectin domain family 4 member L

