



Human CD48 recombinant protein Fc-Tag

OVERVIEW

PRODUCT REFERENCE	PX-P6024
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Ser220
TAG	C-terminal Fc Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	51.57 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P09326
ALIAS/SYNONYMS	Signaling lymphocytic activation molecule 2, TCT.1, CD48 antigen, B-lymphocyte activation marker BLAST-1, BCM1, CD48, BLAST1, Leukocyte antigen MEM-102, SLAMF2, BCM1 surface antigen, SLAM family member 2

