



## Human CD244 recombinant protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P6026
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	Cys22-Pro229
<b>TAG</b>	N-terminal His Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	Escherichia coli (E.coli)
<b>MOLECULAR WEIGHT</b>	25.54 kDa
<b>FORM</b>	
<b>BUFFER</b>	PBS, pH7.5
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	Q9BZW8
<b>ALIAS/SYNONYMS</b>	SLAM family member 4, SLAMF4, h2B4, Signaling lymphocytic activation molecule 4, Natural killer cell receptor 2B4, NKR2B4, NAIL, NK cell type I receptor protein 2B4, CD244, NK cell activation-inducing ligand, 2B4

