



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Human MUC17 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6035
ORIGIN SPECIES	Human
FRAGMENT	Ala4069-Gly4393
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	38.21 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q685J3
ALIAS/SYNONYMS	MUC-3, MUC-17, MUC17, Small intestinal mucin-3, Mucin-17, MUC3



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.