



## Human CD123 / IL3RA recombinant protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P6036
<b>ORIGIN SPECIES</b>	Human
<b>FRAGMENT</b>	Met1-Arg305
<b>TAG</b>	C-terminal His Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Eukaryotic expression
<b>HOST</b>	Mammalian cells
<b>MOLECULAR WEIGHT</b>	36.10 kDa
<b>FORM</b>	
<b>BUFFER</b>	PBS, pH7.5
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	P26951
<b>ALIAS/SYNONYMS</b>	IL-3R subunit alpha, IL-3 receptor subunit alpha, CD123, Interleukin-3 receptor subunit alpha, IL3R, IL-3R-alpha, IL3RA, IL-3RA, Drug Target

