



Human CD213a2 / IL13RA2 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6041
ORIGIN SPECIES	Human
FRAGMENT	Pro238-Glu334
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	13.61 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q14627
ALIAS/SYNONYMS	IL13R, IL-13RA2, IL-13 receptor subunit alpha-2, IL-13R-alpha-2, Interleukin-13-binding protein, CD213a2, IL-13R subunit alpha-2, IL13RA2, Interleukin-13 receptor subunit alpha-2

