



Human CD213a1 / IL13RA1 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6042
ORIGIN SPECIES	Human
FRAGMENT	Gln203-Pro376
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	22.53 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P78552
ALIAS/SYNONYMS	IL13R, IL-13R subunit alpha-1, IL13RA1, IL13RA, IL-13RA1, CD213a1, Interleukin-13 receptor subunit alpha-1, IL-13 receptor subunit alpha-1, CT19, IL-13R-alpha-1, Cancer/testis antigen 19

