



Human CD135 / FLT3 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6043
ORIGIN SPECIES	Human
FRAGMENT	Asn27-Tyr376
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	41.77 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P36888
ALIAS/SYNONYMS	FLT3, Stem cell tyrosine kinase 1, Fetal liver kinase-2, STK-1, STK1, FLK2, CD135, Receptor-type tyrosine-protein kinase FLT3, Fms-like tyrosine kinase 3, FL cytokine receptor, FLK-2, FLT-3

