



## Human OLR1 recombinant protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P6045
<b>ORIGIN SPECIES</b>	Human
<b>FRAGMENT</b>	Leu81-Gln273
<b>TAG</b>	N-terminal His Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Prokaryotic expression
<b>HOST</b>	Escherichia coli (E.coli)
<b>MOLECULAR WEIGHT</b>	24.38 kDa
<b>FORM</b>	
<b>BUFFER</b>	PBS, pH7.5
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	P78380
<b>ALIAS/SYNONYMS</b>	C-type lectin domain family 8 member A, LOX-1, Ox-LDL receptor 1, Lectin-like oxidized LDL receptor 1, OLR1, Oxidized low-density lipoprotein receptor 1, Lectin-like oxLDL receptor 1, LOX1, CLEC8A, hLOX-1, Lectin-type oxidized LDL receptor 1

