



Human CD303 / CLEC4C recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6046
ORIGIN SPECIES	Human
FRAGMENT	Ser67-Ile213
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	19.46 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q8WTT0
ALIAS/SYNONYMS	CLEC4C, C-type lectin superfamily member 7, CLECSF11, C-type lectin domain family 4 member C, Dendritic lectin, CD303, DLEC, HECL, Blood dendritic cell antigen 2, BDCA-2, BDCA2, CLECSF7

