



Human APOC3 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6050
ORIGIN SPECIES	Human
FRAGMENT	Ser21-Ala99
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	21.30 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P02656
ALIAS/SYNONYMS	APOC3, Apolipoprotein C3, Apo-CIII, Apolipoprotein C-III, ApoC-III

