



Human DPEP3 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6052
ORIGIN SPECIES	Human
FRAGMENT	Glu196-Cys488
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	34.69 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9H4B8
ALIAS/SYNONYMS	DPEP3, Dipeptidase 3

