



Recombinant Monkeypox virus/MPXV M1R Protein

OVERVIEW

| | |
|---------------------------|----------------------------------------------|
| PRODUCT REFERENCE | PX-P6056 |
| ORIGIN SPECIES | Monkeypox virus (strain Zaire-96-I-16) (MPX) |
| FRAGMENT | Gly2-Gly183 |
| TAG | C-terminal His Tag |
| DELIVERY CONDITION | Dry Ice |

PRODUCT INFORMATION

| | |
|--------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
| EXPRESSION SYSTEM | Eukaryotic expression |
| HOST | Mammalian cells |
| MOLECULAR WEIGHT | 22.35 kDa |
| FORM | |
| BUFFER | PBS, pH7.5 |
| STABILITY & STORAGE | 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection) |
| UNIPROT | |
| ALIAS/SYNONYMS | M1R, Monkeypox virus/MPXV |

