



Recombinant Monkeypox virus/MPXV HA/B2R Protein

OVERVIEW

PRODUCT REFERENCE	PX-P6062
ORIGIN SPECIES	Monkeypox virus
FRAGMENT	Ser19-Asp274
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	30.83 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	
ALIAS/SYNONYMS	Bifunctional hemagglutinin/type-I membrane glycoproteinImported, Hemagglutinin, MPXV-COP-155, MPXV-SL-155, MPXV-WRAIR155, MPXVgp163, HA, B2R, Monkeypox virus (MPXV)

