



Mouse CTLA-4 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6063
ORIGIN SPECIES	Mouse
FRAGMENT	Ala37-Asp161
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	13.75 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P09793
ALIAS/SYNONYMS	Cd152

