



Mouse PD-1 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6064
ORIGIN SPECIES	Mouse
FRAGMENT	Leu25-Gln167
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	15.73 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q02242
ALIAS/SYNONYMS	Pd1, CD279

