



Human DNMT3A recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6065
ORIGIN SPECIES	Human
FRAGMENT	Ala628-Val912
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	34.98 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	Q9Y6K1
ALIAS/SYNONYMS	Dnmt3a, DNA MTase HsallIA, DNMT3A, M.HsallIA, DNA (cytosine-5)-methyltransferase 3A, DNA methyltransferase HsallIA

