



Human IL-3 recombinant protein

OVERVIEW

| | |
|---------------------------|--------------------|
| PRODUCT REFERENCE | PX-P6067 |
| ORIGIN SPECIES | Human |
| FRAGMENT | Ala20-Phe152 |
| TAG | C-terminal His Tag |
| DELIVERY CONDITION | Dry Ice |

PRODUCT INFORMATION

| | |
|--------------------------------|--|
| EXPRESSION SYSTEM | Prokaryotic expression |
| HOST | Escherichia coli (E.coli) |
| MOLECULAR WEIGHT | 16.71 kDa |
| FORM | |
| BUFFER | PBS, pH7.5 |
| STABILITY & STORAGE | 4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection) |
| UNIPROT | P08700 |
| ALIAS/SYNONYMS | Hematopoietic Growth Factor proteins, IL3, Multipotential colony-stimulating factor, MCGF, Mast cell Growth Factor proteins, Interleukin-3, P-cell-stimulating factor, interleukin 3 |

