



Human CXCL2 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6068
ORIGIN SPECIES	Human
FRAGMENT	Gly34-Asn107
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	10.25 kDa
FORM	
BUFFER	PBS, pH7.5
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P19875
ALIAS/SYNONYMS	GROB, Growth-regulated protein beta, Macrophage inflammatory protein 2-alpha, MIP2A, GRO2, SB-251353, SCYB2, HSF, Gro-beta, Hematopoietic synergistic factor, MIP2-alpha, C-X-C motif chemokine 2, GRO-beta-T

