



Monkeypox virus/MPXV B21R Recombinant Protein, N-His

OVERVIEW

PRODUCT REFERENCE	PX-P6076
ORIGIN SPECIES	Monkeypox virus (strain Zaire-96-I-16) (MPX)
FRAGMENT	Thr291-Thr487
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	24.31 kDa
FORM	
BUFFER	Supplied as solution form in PBS pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)

UNIPROT

ALIAS/SYNONYMS	B21R, Putative membrane-associated glycoprotein, Surface glycoprotein, MPXV-CAM1990_02-171, MPXV-Congo_8-172, MPXV-GAB1988_001-171, MPXV-Ikubi-171, MPXV_RCG2003_358_197, MPXV_SUD2005_01_192, MPXV_ZAI1979_005_197, MPXVgp182, Monkeypox virus (MPXV)
-----------------------	--

