



Monkeypox virus/MPXV C15L Recombinant Protein, C-His

OVERVIEW

PRODUCT REFERENCE	PX-P6079
ORIGIN SPECIES	Monkeypox virus (strain Zaire-96-I-16) (MPX)
FRAGMENT	Met1-Pro176
TAG	C-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	20.46 kDa
FORM	
BUFFER	Supplied as solution form in PBS pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	
ALIAS/SYNONYMS	MPXVgp041, C15L, Monkeypox virus (MPXV)

