



Monkeypox virus/MPXV D6L Recombinant Protein, N-His

OVERVIEW

PRODUCT REFERENCE	PX-P6080
ORIGIN SPECIES	Monkeypox virus (strain Zaire-96-I-16) (MPX)
FRAGMENT	Met1-Lys126
TAG	N-terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Prokaryotic expression
HOST	Escherichia coli (E.coli)
MOLECULAR WEIGHT	16.44 kDa
FORM	
BUFFER	Supplied as solution form in PBS pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	
ALIAS/SYNONYMS	Interleukin-18-binding protein, D6L, Monkeypox virus (MPXV)

