



## Human IL13RA2 recombinant protein

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PX-P6084
<b>ORIGIN SPECIES</b>	Homo sapiens (Human)
<b>FRAGMENT</b>	Met1-Leu342
<b>TAG</b>	C-terminal Fc Tag
<b>DELIVERY CONDITION</b>	Dry Ice

### PRODUCT INFORMATION

---

<b>EXPRESSION SYSTEM</b>	Eukaryotic expression
<b>HOST</b>	Mammalian cells
<b>MOLECULAR WEIGHT</b>	66.05 kDa
<b>FORM</b>	
<b>BUFFER</b>	PBS, pH7.5
<b>STABILITY &amp; STORAGE</b>	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
<b>UNIPROT</b>	Q14627
<b>ALIAS/SYNONYMS</b>	IL13R, IL-13RA2, IL-13 receptor subunit alpha-2, IL-13R-alpha-2, Interleukin-13-binding protein, CD213a2, IL-13R subunit alpha-2, IL13RA2, Interleukin-13 receptor subunit alpha-2

