



Human CD282-TLR2 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6087
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Leu590
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	68.21 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	O60603
ALIAS/SYNONYMS	TIL4, Toll-like receptor 2, TLR2, Toll/interleukin-1 receptor-like protein 4, CD282, Drug Target

