



Human CD61-ITGB3 recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6092
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Pro7-Asp718
TAG	C-Terminal Strep Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	80.12 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P05106
ALIAS/SYNONYMS	GPIIIa, Integrin beta-3, CD61, ITGB3, Platelet membrane glycoprotein IIIa, GP3A, Drug Target

