



Human CD228-MELTF recombinant protein

OVERVIEW

PRODUCT REFERENCE	PX-P6093
ORIGIN SPECIES	Homo sapiens (Human)
FRAGMENT	Met1-Leu738
TAG	C-Terminal His Tag
DELIVERY CONDITION	Dry Ice

PRODUCT INFORMATION

EXPRESSION SYSTEM	Eukaryotic expression
HOST	Mammalian cells
MOLECULAR WEIGHT	79.17 kDa
FORM	
BUFFER	0.01M PBS, pH 7.4.
STABILITY & STORAGE	4°C for short term (1 week), -20°C or -80°C for long term (avoid freezing/thawing cycles; addition of 20-40% glycerol improves cryoprotection)
UNIPROT	P08582
ALIAS/SYNONYMS	CD228, Melanotransferrin, Melanoma-associated antigen p97, MAP97, MF12, MELTF

